

CONTENTS

A.1	MEDIA	148
A.1.1	Artificial Sea Water (ASW)	148
A.1.2	Fe-solution	148
A.1.3	PII metal solution	148
A.1.4	PES-enriched seawater medium (1/3 strength) (Provasoli, 1968)	149
A.1.5	PES-N (ES medium lacking nitrogen)	149
A.2	GENERAL SOLUTIONS	149
A.2.1	General stock solutions	149
A.2.2	Solutions for protein extractions	151
A.2.3	Solutions for the modified Bradford assay for determining protein concentration	151
A.2.4	Solutions for SDS-PAGE	152
A.2.5	Solutions for Coomassie staining of polyacrylamide gels	153
A.2.6	Solutions for 2-DE	154
A.2.7	Western Blot solutions	154

APPENDIX A

MEDIA AND SOLUTIONS

All media were autoclaved at 121 °C for 20 min prior to use, unless otherwise specified.

Water used for making solutions, media and diluting buffers was purified using a Milli-RO Plus (Millipore) water purification system. Ultrapure water used was obtained by further purification of the above water using a Milli-Q Plus (Millipore) water purification system.

A.1 MEDIA

A.1.1 Artificial Sea Water (ASW)

NaCl (Saarchem)	24.7 g
MgCl ₂ .6H ₂ O (Saarchem)	4.7 g
KCl (Saarchem)	0.66 g
CaCl ₂ .2H ₂ O	1.9 g
MgSO ₄ .7H ₂ O (Saarchem)	6.3 g
NaHCO ₃	0.18 g
Water to	1 L

Autoclave

A.1.2 Fe-solution

Fe(NH ₄) ₂ (SO ₄) ₂ .6H ₂ O	702 mg
Na ₂ EDTA	600 mg
Water to	1 L

A.1.3 PII metal solution

Na ₂ EDTA	100 mg
H ₃ BO ₃	114 mg
FeCl ₃ .6H ₂ O	4.9 mg
MnSO ₄	16.4 mg
ZnSO ₄ .7H ₂ O	2.2 mg
CoSO ₄ .7H ₂ O	0.48 mg

Water to	100 ml
----------	--------

A.1.4 PES-enriched seawater medium (1/3 strength) (Provasoli, 1968)

NaNO ₃	350 mg
Na ₂ glycerophosphate 5H ₂ O	50 mg
Fe solution	25 ml
PII	25 ml
Vitamin B ₁₂	10 µg
Thiamine	0.5 mg
Biotin	5 µg
Tris buffer (Sigma Co.)	500 mg
Water to	100 ml

Adjust pH to 7.8, autoclave and store at 10 °C.
Add 6.6 ml to 1 L ASW.

A.1.5 PES-N (ES medium lacking nitrogen)

Prepared exactly the same as ES medium except no NaNO₃ is added and Fe₂(SO₄) instead of Fe(NH₄)₂(SO₄)₂·6H₂O when making Fe-solution.

A.2 GENERAL SOLUTIONS

A.2.1 General stock solutions

- 1 N NaOH**

NaOH (Saarchem)	4 g
water to	100 ml

Store in a plastic bottle.

- 1 M HCl**

37% HCl (Saarchem)	8.4 ml
water to	100 ml

Store in a foil covered glass bottle.

- 0.5 M EDTA**

Na ₂ EDTA (Saarchem)	186.1 g
water to	800 ml

The pH was adjusted to 8.0 with NaOH pellets prior to making the volume to 1 l with water and autoclaving.

- **1 M Tris-HCl**

Tris base	121.5 g
water	800 ml

The pH was adjusted to 8.0 with HCl prior to making the volume to 1 l with water and autoclaving.

- **25% (w/v) Sodium dodecyl sulphate (SDS)**

SDS	50 g
sterile water	200 ml

Stir on warm plate to dissolve. Do not autoclave. The concentration of the SDS solution can be altered by altering the amount of SDS added.

- **70% (v/v) Ethanol (EtOH)**

Absolute ethanol	70 ml
water to	100 ml

Do not autoclave.

- **50% (v/v) Glycerol**

glycerol	25 ml
water to	50 ml

- **2 M Sodium phosphate, mono-sodium (NaH_2PO_4)**

$\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$	55.2 g
water to	200 ml

- **2 M Sodium phosphate, di-sodium (Na_2HPO_4)**

$\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$	71.2 g
water to	200 ml

- **1 M Phosphate buffer (pH 6.8)**

$\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ (2 M)	39 ml
$\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$ (2 M)	61 ml
water to	200 ml

- **0.1 M Phosphate buffer**

Phosphate buffer (1 M, pH 6.8)	10 ml
Water to	100 ml

- **DTT (30%)**
0.3 g DTT added to 1 ml distilled water
50 µl aliquots were made and stored at -20°C
- **ASB 14 (10%)**
0.1 g ASB 14 added to 1 ml distilled water
50 µl aliquots were made and stored at -20°C

A.2.2 Solutions for protein extractions

- **Protein extraction buffer**

Tris (0.5 M) pH 7.5	5 ml
EDTA (0.5 m)	200 µl
Triton	1%
β-mercaptoethanol	2%
water to	10 ml

- **0.1 M ammonium acetate in methanol (w/v)**

3.85 g ammonium acetate added to 500 ml of methanol

- **80% acetone (v/v)**

80ml of acetone added to 20ml of water

- **Urea Lysis Buffer (ULB)**

Urea (8 M)	24.0 g
------------	--------

Dissolve in 25 ml water by stirring and then add:

CHAPS (2 %) (Calbiochem)	0.25 g
Thiourea (Sigma-Aldrich)	7.6 g

Make up to 50 ml with distilled water and filter sterilize.

A.2.3 Solutions for the modified Bradford assay for determining protein concentration

- **Bovine serum albumin (BSA) (1 mg/ml⁻¹)**

BSA (Roche)	0.01 g
ULB to	10 ml

- **Bradford's reagent**

Dye Reagent Concentrate (Bio-Rad)	1 ml
sterile water to	5 ml

Make up fresh before every use

- **0.1 M HCl**

HCl	10 ml
Water to	100 ml

A.2.4 Solutions for SDS-PAGE

- **Stacking gel**

4X Tris-Cl / SDS solution (pH 6.8)

Tris base	6.05 g
SDS	0.4 g

Adjust pH to 6.8 with HCl
 Add water to 100 ml
 Filter-sterilize and store at room temperature

Separating gel

4X Tris-Cl / SDS solution (pH 8.8)

Tris base (1.25 M)	91.0 g
SDS	2.0 g

Adjust pH to 8.8 with 1 M HCl
 Add water to 500 ml
 Filter-sterilize and store at room temperature

- **10 x SDS-PAGE running buffer**

Glycine	150 g
Tris	30 g
SDS	20 g
water to	1000 ml

- **1x SDS PAGE running buffer**

Dilute 10 x SDS-PAGE running buffer 1:10 with water before use

- **10% (w/v) Ammonium persulfate (AMPS)**

AMPS	1 g
Water to	10 ml

Do not autoclave. Filter sterilize before aliquoting into sterile microfuge tubes and store at -20°C.

- **5x SDS-PAGE Sample Application Buffer (SAB)**

Tris-Cl (pH 6.8)	250 mM
DTT (Promega)	500 mM
SDS (Sigma)	10%
Bromophenol blue (Saarchem)	0.5%
Glycerol	10 ml
water to	5 ml

A.2.5 Solutions for Coomassie staining of polyacrylamide gels

- **Coomassie SDS-PAGE gel stain**

Coomassie brilliant blue R250 (Sigma)	1 g
Methanol (50%)	225 ml
Glacial acetic acid (10%)	50 ml
water to	

Dissolve the Coomassie Brilliant Blue R in methanol before adding acetic acid and water

- **Coomassie SDS-PAGE destain**

Methanol	50 ml
Glacial acetic acid	70 ml
water to	1000 ml

- **SDS-PAGE Gel storage solution**

Glacial acetic acid	7 ml
Water to	100 ml

A.2.6 Solutions for 2-DE

- **Rehydration solution**

ULB	X ml
DTT	2%
ASB 14	0.4%
Biolyte 3-10 ampholytes (20%)	1%
Trace amount of bromophenol blue	
Protein sample (250 µg)	Y µl

Made in a total volume of 140 µl.

- **Equilibration buffer I (EQB I)**

Urea (6M)	72.1 g
Tris-HCl (pH 8.8, 0.375 M)	50 ml
Glycerol	20%
SDS	2%
Bromophenol blue	0.0005%

Made up to 200 ml with distilled water, aliquoted into 10 ml volumes and stored at -20°C.

0.1 g DTT was added to 10 ml of EQB I prior to the equilibration step II.

- **Equilibration buffer II (EQB II)**

Made up as for EQB I except that 0.48g of iodoacetamide was added to EQB II prior to equilibration step I.

- **0.5% Agarose**

0.5 g agarose in 1 X SDS running buffer with a trace amount of bromophenol blue.

A.2.7 Western Blot solutions

- **Ponceau S**

Ponceau S powder	0.1%
Ethanol	5%
Made in distilled water	

This stain is re-usable but due to its sensitivity to light needs to be covered in foil and stored and used in the dark.

- **10X TBS Stock**

50 mM Tris	30.25g
150 Mm NaCl	43.8g

Distilled water to 300 ml
 Adjust pH to 7.4 with HCl
 Distilled water to 500ml
 Autoclave

- **Blocking buffer (5 % skim milk)**

Skim-milk (Elite)	5 g
10 X TBS (1 X)	10 ml

Water to	100 ml
----------	--------

- **TBST (0.1 % Tween in 1X TBS)**

Tween 20 (Saarchem)	0.1 ml
10X TBS (1 X)	100 ml
Water to	1000 ml

Make up fresh
 Do not autoclave

- **Towbin Buffer**

Tris base	3.03 g
Glycine (Merck)	14.42 g
Methanol (20%) (Merck)	200 ml
Water to	1000 ml